


**NANJING PERLOVE MEDICAL EQUIPMENT CO.,LTD.**

Add: No.168 Kaiyuan Road, Jiangning District, Nanjing, 211100, China  
Tel: +86-25-87187780 +86-25-68571666-9350 9353  
Fax: +86-25-87187780  
Http://www.perlove.com.cn  
E-mail: overseas@perlove.com.cn



**OUR UPCOMING EXHIBITIONS:**

You can see more details on our website  
Http://www.perlove.com.cn



**PERLOVE MEDICAL**

**For Overseas Market:** Nanjing Perlove Medical Equipment Co.,Ltd.

**For Domestic Market:** Nanjing Perlove Medical Equipment Co.,Ltd.

Beijing Per-Builders Science and Technology Co.,Ltd.

Zuhai Perlead Medical Equipment Co.,Ltd.

# VET1100

## Vet Digital Radiography System

Version 1.0 Year 2017

Http://www.perlove.com.cn



# Veterinary Digital Radiography System

## Advanced Technology

### Innovative

- ☺ Mature flat-panel detector ensures excellent image quality and reliable clinical diagnosis.
- ☺ Fast imaging time, high efficiency, directly imaging in computer without film processing.



## A-Si Flat Panel Detector

- Pixel pitch 154um
- Effective area 17\*17inch
- Spatial resolution 3.6Lp/mm
- A/D 14bit
- Model Amorphous Silicon
- Pixel matrix 3072\*3072



Size 17\*17M

## Optimized Human-friendly Design

Multiple self-protection and fault alarming function, safe and reliable.

Automatically save parameters when suddenly power off, no loss of data.

Four way floating table, electromagnetic brake ensures accurate positioning.

## Simple Operation Sharp Image

Special veterinary radiography modes, easy to operate and safe use.

DICOM 3.0 interface supports concise image processing and work list output.

Configuration VET1100

Clinical application Radiography for limb abdomen chest etc

Output power 4.5KW

Inverter frequency 40kHz

Tube Voltage 40kV-120kV

HF High Voltage Generator and Tube Tube Current 20mA-100mA

Radiography(mAs) 1.0mAs-180mAs

Power supply 220V

Exposure method Line control and remote control



## Unique Clinical Experience

High quality X-ray monoblock assures more details, less dose.

Medical fiberboard table surface, scratch-prevented.